COTININ TEST (COT)
008A086

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Acetaminophen	Diazepam	(+/-)-Epinephrine	Penicillin-G	Tyramine
Amitriptyline	Dopamine	Chlorprothixene	Pentobarbital	Vitamin C
D-Amphetamine	Doxylamine	r-Cyclodextrin	Pheniramine	Theophylline
Ampicillin	Ecgonine methyl ester	Cyclobenzaprine	Phenothiazine	p-Naphthaleneacetic acid
Aspartame	(-)-Ephedrine	5,5-Diphenylhydantoin	L-Phenylephrine	D(+)-Trehalose
Aspirin	(+/-)-Ephedrine	Nicotinic	4-Dimethylaminoantipyrine	(+/-)-Soprotrenol
Atropine	Erythromycin	Methamphetamine	p-Phenylethylamine	(+)-Ephedrine
Benzocaine	EDDP	Methaqualone	Phentermine	(+/-)-Phenylpropanolamine
Benzoyllecgonine	Furosemide	Methadone	Procaine	Methapyrilene
[+]-Brompheniramine	Glucose	(+/-)-3,4-MDMA	Promethazine	(+/-)-Phenylpropanolamine
Caffeine	Guaiaacol glyceryl ether	Morphine	d-Propoxyphene	Methapyrilene
Chloroquine	Hydrocodone	Morphine-3+d-glucuronide	Secobarbital	(+/-)-Phenylpropanolamine
(+)-Chlorpheniramine	Hydromaphone	(1R,2S)-(-)-N-Methyl-Ephedrine	Sodium Chloride	Acid Acetylsalicylic Acid
(+/-)-Chlorpheniramine	Lidocaine	Naloxone	Tenocyclidine	Triprolidine Hydrochloride
Codeine	Meperidine	Naltrexone	11-nor- Δ^9 -Tetrahydrocannabinol	
			- β -carboxylic acid	
Creatine	Methadol	(+)-Naproxen	Thioridazine	
(-)-Deoxyephedrine	Maprotiline	Nortriptyline	Trimethobenzamide	
Dextromethorphan	DL-Homatropine	Oxalic Acid	Trifluoperazine	

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- Zevin S, Jacob P and Benowitz N. Cotinine effects on nicotine metabolism. *Clinical Pharmacology & Therapeutics*, 1997 June 61(6): 649-54.
- Zuccaro P, Pichini S, Altieri I, Rosa M, Pellegrini M and Pacifici R. Interference of nicotine metabolites in cotinine determination by RIA. *Clinical Chemistry*, 1997 Jan 43(1): 180-1.
- Benowitz NL. Cotinine as biomarker of environmental tobacco smoke exposure. *Epidemiologic Reviews*, 1996 18(2): 188-204.
- Curvall M and Vala FK. Nicotine and metabolites: analysis and levels in body fluids. Nicotine and related alkaloids: absorption, distribution, metabolism and excretion. Edited by Gorrod JW and Wahren J. Published in 1993 by Chapman & Hall, London. ISBN 0 142 55740 1.

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Distributor pro ČR a SR: JK-Trading. spol. s r.o., Za Opravnou 5, 150 00 Praha 5, tel: 257 220 760, email: praha@jktrading.cz, www.jktrading.cz

ulti med Products (Deutschland) GmbH
Reeshoop 1 • 22926 Ahrensburg
Telefon: 04102 - 80090
Fax: 04102 - 50082
e-mail: info@ultimed.de

ulti med Products (Belgium) BVBA
Honzebroekstraat 137 • 8800 Roeselare
Phone: +32 +51 200 425
Fax: +32 +51 200 449
e-mail: belgium@ultimed.com