NADAL® Syphilis Test (test cassette)

REF 203902N-20

| | | | | | |
|-----------|----------------------------|----|-----------|--------------------------|----|
| de | Gebrauchsanweisung | 2 | pt | Instruções de Utilização | 26 |
| en | Instructions for use | 6 | cs | Návod k použití | 30 |
| fr | Instructions d'utilisation | 10 | | Symbols | 35 |
| es | Instrucciones de uso | 14 | | Our Teams | 36 |
| it | Istruzioni per l'uso | 18 | | | |
| pl | Sposób użycia | 22 | | | |

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Carl-Zeiss-Strasse 12
47445 Moers
Germany

Moers
Tel: +49 (2841) 99820-0
Fax: +49 (2841) 99820-1

Regensburg
Tel: +49 941 29010-0
Fax: +49 941 29010-50

www.nal-vonminden.com
info@nal-vonminden.com

Directors:
Sandra von Minden
Roland Meißner
Thomas Zander

Commercial reg. Kleve
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

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15. Referenzen

- Centers for Disease Control (CDC). Chlamydia trachomatis infections. Policy guidelines for prevention and control. MMWR Morb Mortal Wkly Rep. 1985 Aug 23; 34 Suppl 3: 535-745.
- Tichonova L, Borisenko K, Ward H, Meheus A, Gromyko A, Renton A. Epidemics of syphilis in the Russian Federation: trends, origins, and priorities for control. Lancet. 1997 Jul 19; 350(9072): 210-3.
- Norgard MV, Chamberlain NR, Swancutt MA, Goldberg MS. Cloning and expression of the major 47-kilodalton surface immunogen of *Treponema pallidum* in *Escherichia coli*. Infect Immun. 1986 Nov; 54(2): 500-6.

1. Intended Use

The NADAL® Syphilis Test is a lateral flow chromatographic immunoassay for the qualitative detection of anti-*Treponema pallidum* (*T. pallidum*) IgM and IgG in human whole blood, serum or plasma specimens. The test is intended for use as an aid in the diagnosis of syphilis in patients suspected of having a syphilis infection (see section 12 'Limitations'). The test procedure is not automated and requires no special training or qualification. The NADAL® Syphilis Test is designed for professional use only.

2. Introduction and Clinical Significance

Treponema pallidum (*T. pallidum*), a spirochete bacterium with an outer cell membrane and a cytoplasmic membrane, is the causative agent of the venereal disease known as syphilis. Although syphilis rates are declining in the United States after an epidemic between 1986 and 1990, the incidence of syphilis in Europe has been increasing since 1992, especially in the states of the Russian Federation, where peaks of 263 cases per 100,000 have been reported. In addition, the positive rate of serological test results for syphilis in HIV-infected individuals has recently been on the rise.

The serological detection of specific antibodies to *T. pallidum* has been long recognised in the diagnosis of syphilis since the natural course of the infection is characterised by periods without clinical manifestations. The antibody response to *T. pallidum* can be detected within 4 to 7 days after the syphilis chancre appears, allowing the early detection and diagnosis of a syphilis infection.

A variety of antigens, such as Cardiolipin (RPR test), VDRL antigen and *T. pallidum* extracts derived from *in-vitro* culture or inoculated rabbit testes, have been used in syphilis serological tests. However, RPR and VDRL antigens are not treponemal specific and whole *T. pallidum* extracts are not reproducible and contain a certain amount of contaminating materials such as flagella, which may lead to a non-specific reaction in assays of test serum.

3. Test Principle

The NADAL® Syphilis Test enables the detection of anti-*T. pallidum* IgM and IgG through the visual interpretation of colour development on the internal test strip. Specific recombinant *T. pallidum* antigens are immobilised in the test line region (T) of the membrane. During the test, the specimen reacts with the *T. pallidum* antigens which are conjugated to coloured particles and precoated onto the conjugate pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there are a sufficient number of anti-*T. pallidum* antibodies in the specimen, a coloured line will develop in the test line region (T) of the membrane. The presence of this coloured line indicates a positive result, while its absence indicates a negative result. The formation of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® Syphilis test cassettes, incl. disposable pipettes (40 µL)
- 1 buffer (4 mL)*
- 1 package insert

*containing the following preservative: ProClin™ 300: <0.03%.

No hazard labelling for ProClin™ 300 is required according to Regulation (EC) № 1272/2008 CLP. Concentrations are below the exemption threshold of <0.03%.

5. Additional Materials Required

- Specimen collection containers (suitable for specimen material to be tested)
- Centrifuge (for serum or plasma specimens only)
- Alcohol pads
- Lancets (for fingerstick whole blood specimens only)
- Timer

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the entire instructions for use prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.

- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

8. Specimen Collection and Preparation

The NADAL® Syphilis Test can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.

To collect fingerstick whole blood specimens:

- Wash the patient's hand with soap and warm water or clean it with an alcohol pad. Allow it to dry.
- Massage the hand, without touching the puncture site, by rubbing along the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first drop of blood.
- Gently rub the hand from the wrist to the palm, and then to the finger to form a rounded drop of blood over the puncture site.

Fingerstick whole blood should be tested immediately.

Venipuncture whole blood specimens

Containers containing anticoagulants, such as EDTA, sodium citrate or sodium heparin should be used for the preparation of venous whole blood or plasma specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time.

If the test is to be run within 2 days of specimen collection, whole blood collected by venipuncture should be stored at 2-8°C.

Do not freeze whole blood specimens.

Serum and plasma specimens

Separate serum or plasma from blood as soon as possible to avoid haemolysis. Use only clear, non-haemolysed specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum and plasma specimens can be stored at 2-8°C for up to 3 days. For long-term storage, specimens should be kept at -20°C.

Bring specimens to room temperature prior to testing. Frozen specimens should be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.

If specimens are to be shipped, they should be packed in compliance with all applicable regulations for the transportation of etiologic agents.

Icteric, lipemic, haemolysed, viscous, heat-treated and contaminated specimens may lead to inaccurate test results.

9. Test Procedure

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.

2. Place the test cassette on a clean and level surface.

3. a) For serum, plasma or venous whole blood specimens


- Holding a pipette vertically, transfer 2 drops (approximately 80 µL) of the specimen to the specimen well (S) of the test cassette.
- Holding the buffer bottle vertically, add 1 drop (approximately 40 µL) of buffer to the specimen well (S).

Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.

b) For fingerstick whole blood specimens using hanging drops:

- Position the patient's finger so that a drop of blood is exactly above the specimen well (S) of the test cassette.
- Allow 2 hanging drops of fingerstick whole blood to fall into the centre of the specimen well (S) of the test cassette.
- Holding the buffer bottle vertically, add 1 drop (approximately 40 µL) of buffer to the specimen well (S). **Avoid squeezing the finger, as this might lead to inaccurate test results.**

4. Start the timer.

As the test begins to run, you will observe a coloured liquid migrate along the membrane.

- Wait for the coloured line(s) to appear. Read the test result after 10 minutes. Do not interpret the result after more than 30 minutes.

10. Result Interpretation

Positive:

A coloured line develops in the control line region (C) and another coloured line develops in the test line region (T).

Note: The colour intensity in the test line region (T) may vary depending on the analyte concentration present in the specimen. Any shade of colour in the test line region (T) should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

Negative:

A coloured line develops in the control line region (C). No line develops in the test line region (T).

**Invalid:**

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

The internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal positive procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® Syphilis Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection of anti-*T. pallidum* IgM and IgG in human whole blood, serum or plasma specimens only.
- Neither the quantitative value nor the rate of increase/decrease in the concentration of *T. pallidum* antibodies can be determined using this qualitative test.
- The accuracy of the test depends on the quality of the specimen. Inaccurate test results may occur due to improper specimen collection or storage (see section 8 'Specimen Collection and Preparation').
- Some specimens containing unusually high titers of heterophile antibodies or rheumatoid factor may affect test results.
- As with all diagnostic tests, all results should be interpreted by a physician in conjunction with other available clinical information.
- The NADAL® Syphilis Test only detects the presence of *T. pallidum* antibodies in specimens and should not be used as the sole criterion for a diagnosis of *T. pallidum* infection.
- If the test result is negative but clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of a *T. pallidum* infection.

13. Performance Characteristics**Clinical performance****Diagnostic sensitivity and specificity**

The NADAL® Syphilis Test was evaluated using clinical whole blood specimens in comparison to another commercially available TPHA Syphilis rapid test.

The results are presented in the following table:

| NADAL® Syphilis Test | Another TPHA Syphilis rapid test | | |
|----------------------|----------------------------------|----------|-------|
| | Positive | Negative | Total |
| | Positive | 394 | 4 |
| | Negative | 2 | 540 |
| | Total | 396 | 544 |

Diagnostic sensitivity: 99.5% (98.2% - 99.9%)*

Diagnostic specificity: 99.3% (98.1% - 99.8%)*

Overall agreement: 99.4% (98.6% - 99.8%)*

*95% confidence interval

Analytical performance**Analytical specificity****Cross-reactivity study**

Negative whole blood specimens were spiked with the following, potentially cross-reacting substances and tested using the NADAL® Syphilis Test:

| Substance | Concentration |
|------------------------------------|---------------|
| HAMA | 1000 ng/mL |
| HBsAb | Positive |
| HBcAb | Positive |
| anti- <i>H. pylori</i> IgG | Positive |
| anti-rubella virus IgM | Positive |
| Rheumatoid factors (RF) | 8400 IU/mL |
| anti-HBcAg IgM | Positive |
| anti-HCV IgM | Positive |
| anti-EBV Ab | Positive |
| anti- <i>Toxoplasma gondii</i> IgM | Positive |
| HBsAg | Positive |
| HBeAb | Positive |
| anti-HIV IgM | Positive |
| anti-CMV IgM | Positive |

No cross-reactivity with the specimens was observed when tested using the NADAL® Syphilis Test.

Interference study

Negative whole blood and low positive serum specimens spiked with the following potentially interfering substances at the concentrations listed below showed no interference with the NADAL® Syphilis Test.

| Substance | Concen- tra-tion | Substance | Concentration |
|----------------------|---------------------|------------------|---------------|
| Acetaminophen | 0.2 mg/mL | Caffeine | 0.2 mg/mL |
| Acetylsalicylic acid | 0.2 mg/mL | Gentisic acid | 0.2 mg/mL |
| Ascorbic acid | 20 mg/mL | Albumin | 20 mg/mL |
| Creatine | 2 mg/mL | Haemo- globin | 0.011mg/mL |
| Bilirubin | 10 mg/mL | Oxalic acid | 6 mg/mL |

Precision

Repeatability

Repeatability was established by testing 20 replicates of negative, low, medium and high positive specimens using 3 lots of the NADAL® Syphilis tests. >99% of the specimens were correctly identified (20/20 correct tests per concentration, 95% confidence interval: 95.5-100%). The NADAL® Syphilis Test demonstrated acceptable repeatability.

Reproducibility

Reproducibility was established by testing 5 replicates of negative, low, medium and high positive specimens. Testing was performed by 4 operators using 3 independent NADAL® Syphilis test lots at 3 sites on 5 separate days. >99% of the specimens were correctly identified (100/100 correct tests per concentration, 95% confidence interval: 99.3%-100%). The NADAL® Syphilis Test demonstrated acceptable reproducibility.

14. Serious incident reporting

In the case of any serious incidents related to the performance of the NADAL® Syphilis Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

15. References

1. Centers for Disease Control (CDC). Chlamydia trachomatis infections. Policy guidelines for prevention and control. MMWR Morb Mortal Wkly Rep. 1985 Aug 23; 34 Suppl 3: 535-745.
2. Tichonova L, Borisenko K, Ward H, Meheus A, Gromyko A, Renton A. Epidemics of syphilis in the Russian Federation: trends, origins, and priorities for control. Lancet. 1997 Jul 19; 350(9072): 210-3.
3. Norgard MV, Chamberlain NR, Swancutt MA, Goldberg MS. Cloning and expression of the major 47-kilodalton surface immunogen of *Treponema pallidum* in *Escherichia coli*. Infect Immun. 1986 Nov; 54(2): 500-6.

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15. Bibliografia

1. Centers for Disease Control (CDC). Chlamydia trachomatis infections. Policy guidelines for prevention and control. MMWR Morb Mortal Wkly Rep. 1985 Aug 23; 34 Suppl 3: 53S-74S.
2. Tichonova L, Borisenko K, Ward H, Meheus A, Gromyko A, Renton A. Epidemics of syphilis in the Russian Federation: trends, origins, and priorities for control. Lancet. 1997 Jul 19; 350(9072): 210-3.
3. Norgard MV, Chamberlain NR, Swancutt MA, Goldberg MS. Cloning and expression of the major 47-kilodalton surface immunogen of *Treponema pallidum* in *Escherichia coli*. Infect Immun. 1986 Nov; 54(2): 500-6.

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