NADAL[®] H. pylori Ab Test (test cassette)

REF 262001

DE	Gebrauchsanweisung	2	CZ	Návod k použití	23
EN	Instructions for use	5	NO	Bruksanvisning	26
FR	Instructions d'utilisation	8	NL	Gebruiksaanwijzing	29
ES	Instrucciones de uso	11		Symbols	35
IT	Istruzioni per l'uso	14		Our Teams	36
PL	Sposób użycia	17			
PT	Instruções de Utilização	20			

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Rev. 0, 2017-08-21 OM/KB

1. Intended Use

The NADAL® H. pylori Ab Test is a rapid visual immunoassay for the qualitative presumptive detection of specific IgG against *Helicobacter pylori* (*H. pylori*) in human whole blood, serum or plasma specimens. The test is intended for use as an aid in the diagnosis of *H. pylori* infections. The NADAL® H. pylori Ab Test is designed for professional use only.

2. Introduction and Clinical Significance

Gastritis and peptic ulcers are among the most common human diseases. Since the discovery of the bacterium *H. pylori* (Warren & Marshall, 1983), many reports have suggested that this bacterium is one of the major causes of ulcer diseases (Anderson & Nielsen, 1983; Hunt & Mohamed, 1995; Lambert et al., 1995). Though the exact role of *H. pylori* in the pathogenesis is not yet fully understood, eradication of *H. pylori* from the gastrointestinal tract is associated with the healing of ulcer diseases.

The human serological responses to infections with *H. pylori* have been demonstrated (Varia & Holton, 1989; Evans et al., 1989). The detection of specific anti-*H. pylori* IgG has been demonstrated to be an accurate method of detecting *H. pylori* infections in symptomatic patients. It should be noted, however, that *H. pylori* may also colonize some asymptomatic people. A serological test can therefore be used either as a supplement to endoscopy or as an alternative measure in symptomatic patients.

3. Test Principle

The NADAL® H. pylori Ab Test detects specific anti-*H. pylori* IgG through visual interpretation of colour development on the internal test strip. *H. pylori* antigens are immobilized in the test line region of the membrane. During testing, the specimen reacts with further *H. pylori* antigens conjugated to coloured particles and precoated onto the conjugate pad of the test. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there are sufficient anti-*H. pylori* antibodies in the specimen, a coloured line will appear in the test line region of the membrane. The appearance of a coloured line in the control line region serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® H. pylori Ab test cassettes incl. disposable pipettes
- 1 bottle containing Whole Blood Buffer
- 1 package insert

5. Additional Materials Required

- Specimen collection container
- Timer
- Centrifuge (for serum collection)

6. Storage & Stability

The test kit should be stored at 2-30°C until the expiry date printed on the sealed foil pouch. The test cassette should remain in the sealed foil pouch until use. Do not freeze the test kit. Care should be taken to protect the components of the test kit from contamination. Do not use the test if there is

evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to false results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test after the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Humidity and temperature can adversely affect test results.
- Used testing materials should be discarded according to local regulations.

8. Specimen Collection and Preparation

The NADAL® H. pylori Ab Test is intended for use with human whole blood, serum or plasma specimens only.

Only clear, non-haemolysed specimens are recommended for use with this test.

Separate serum or plasma as soon as possible to avoid haemolysis.

Perform testing immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum and plasma specimens can be stored at 2-8°C for up to 3 days. For long-term storage, specimens should be kept at temperatures below -20°C.

Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.

For whole blood storage specimen collection containers containing anticoagulants such as EDTA, citrate or heparin should be used.

Bring specimens to room temperature prior to testing. Frozen specimens should be completely thawed and well mixed prior to testing. Avoid repeated freeze-thaw cycles.

If specimens are to be shipped, pack them in compliance with all applicable regulations for transportation of etiological agents.

Icteric, lipemic, haemolysed, heat-treated and contaminated specimens may lead to inaccurate test results.

9. Test Procedure

Bring tests, Whole Blood Buffer, specimens and/or controls to room temperature (15-30°C) before use. The test cassette should be used immediately after opening the foil pouch.


1. Remove the test cassette from the sealed foil pouch and place it on a clean and level surface. Label the test cassette with the patient or control identification. For the best results, the assay should be performed within one hour.

2 a) Transfer 3 drops of serum or plasma specimen (approximately 75 µL) into the specimen well (S) of the test cassette with the provided disposable pipette.

OR

b) Transfer 2 drops of whole blood specimen (approximately 50 µL) into the specimen well (S) of the test cassette with the provided disposable pipette, then add 1 drop of Whole Blood Buffer.

OR

c) Allow 2 drops of fingerstick whole blood specimen (approximately 50 µL) to fall into the center of the specimen well (S) of the test cassette, then add 1 drop of Whole Blood Buffer.

Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.

As the test begins to run, you will observe a coloured liquid migrate along the membrane.

3. Start the timer.

4. Wait for the coloured line(s) to appear. Read the test result after 10 minutes. Do not interpret the result after more than 20 minutes.

10. Result Interpretation

Positive

Two coloured lines appear on the membrane. One line appears in the control line region (C) and the other line appears in the test line region (T).

Negative

Only one coloured line appears in the control line region (C). No coloured line appears in the test line region (T).

Invalid

The control line fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Note:

The colour intensity in the test line region (T) may vary depending on the concentration of the analyte present in the specimen. Therefore, any shade of colour in the test line region should be considered positive. Note that this is a qualitative test only and it cannot determine the concentration of the analyte in the specimen.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure.

11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® H. pylori Ab Test is for professional *in-vitro* diagnostic use only and should solely be used for the qualitative detection of anti-*H. pylori* antibodies. No meaning should be inferred from the colour intensity or width of any apparent lines.
- This test should be used for symptomatic patients with gastrointestinal disorders. A diagnosis of gastritis or peptic ulcers should be made in conjunction with other clinical and laboratory findings.
- A positive result suggests only the presence of specific anti-*H. pylori* antibodies and does not distinguish between active and past infections. A positive result is not necessarily indicative of gastrointestinal disease.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time rule out the possibility of *H. pylori* infection, as the concentration of anti-*H. pylori* antibodies may be below the detection level of the test.

13. Performance Characteristics

Clinical Performance

A total of 617 specimens from both symptomatic and asymptomatic individuals were tested using the NADAL® H. pylori Ab Test. A biopsy (bacterial culture), a histology and a rapid urease test (RUT) were used as a reference method. The results are presented in the following table.

NADAL® H.pylori Ab Test – Biopsy/ Histology/ RUT

		NADAL® H. pylori Ab Test		
		+	–	Total
Biopsy/ Histology/ RUT	+	246	18	264
	–	10	343	353
	Total	256	361	617

Relative Sensitivity: 93.2% (89.5%-95.9%)*

Relative Specificity: 97.2% (94.9% - 98.6%)*

Overall agreement: 95.5% (93.5% - 97.0%)*

*95% Confidence Interval

Interference

The interference of known, potentially interfering substances on the performance of the NADAL® H. pylori Ab Test was studied by spiking anti-*H. pylori* IgG negative and positive blood specimens with these substances (see the table below). At the concentrations tested, no interference was observed.

Substance	Concentration
Haemoglobin	2 g/l
Bilirubin	342 µmol/l
Triglycerides	37 mmol/l
Albumin	60 g/l
γ globulin	60 g/l
HAMA	780.73 IU/ml
Rheumatoid factor	1035 IU/ml
Metronidazole	701 µmol/l
Quinine	148 µmol/l
Rifampicin	78.1 µmol/l
Aspirin	4.34 mmol/l
Paracetamol	199 µmol/l
Ibuprofen	2425 µmol/l
Caffeine	308 µmol/l
Ethanol	86.8 mmol/l

Cross-reactivity

Cross-reactivity with other immunoproteins and micro-organisms that might be present in human blood specimens has been studied. The NADAL® H. pylori Ab Test shows no cross-reactivity with specimens which tested positive for

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nal von minden GmbH