NADAL[®] Trichomonas Vaginalis Test (test cassette)

REF 840003N-10

DE	Gebrauchsanweisung	2	PL	Sposób użycia	18
EN	Instructions for use	6	CZ	Návod k použití	22
FR	Instructions d'utilisation	10		Symbols	27
ES	Instrucciones de uso	14		Our Teams	28

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- World Health Organisation (2001) Global prevalence and incidence of selected curable sexually transmitted infections: overview and estimates. World Health Organisation, Geneva, Switzerland.
- Ledru S, Meda N, Fofana M, Soula G, Bazie AJ, Chiron JP (1992) Etiologic study of genitourinary infections in women of child bearing age in Bobo-Dioulasso, Burkina Faso Sex Transm Dis; 23(2):151-156.
- Yereli K, Balcioglu IC, Degerli K, Ozbilgin A, Daldal N (1997) Incidence of Trichomonas vaginalis among women having vaginal discharge, in Manisa, Turkey J Egypt Soc Parasitol 27 (3): 905 – 911.

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1. Intended Use

The NADAL® *Trichomonas Vaginalis* Test is a simple, one step chromatographic immunoassay for the rapid, qualitative detection of *Trichomonas vaginalis* in vaginal swab specimens. This test is intended for use as an aid in the diagnosis of trichomoniasis and designed for professional use only.

2. Introduction and Clinical Significance

Trichomoniasis is a common sexually transmitted infection, which occurs in both men and women, though symptoms are more common in women. It is caused by a protozoan parasite called *Trichomonas vaginalis*. *Trichomonas vaginalis* is pear-shaped and 7-30 µm long. It has five flagella. The large nucleus is usually located at the wider, anterior end and contains many chromatin granules and a small karyosome. The cytoplasm also contains many granules, but these are often not seen in Giemsa-stained specimens.

The most common symptoms of trichomoniasis include:

- Discomfort during intercourse and urination (in both sexes)
- Frothy yellow-green vaginal discharge with a strong, unpleasant odour (in women)
- Irritation and itching of the genital area (in women)
- Rarely, trichomoniasis is associated with lower abdominal pain (in women)

In women, trichomoniasis infection is usually confined to the vagina. However occasionally the microorganisms can also invade the urinary tract and cause cystitis (inflammation of the bladder). Trichomoniasis can be severe enough to cause an abnormal Pap smear, from which the microorganisms also can be detected, and can recur. Although *Trichomonas vaginalis* can cause small red lesions on the cervix, it does not invade the uterus or fallopian tubes, nor does it affect fertility.

There is evidence to suggest that individuals with trichomoniasis are at a higher risk for contracting HIV infection during sexual contact with an HIV-infected partner.

3. Test Principle

The NADAL® *Trichomonas Vaginalis* Test is based on the principle of highly specific binding reactions between antigens and antibodies, which are used for the analysis of specific substances in specimens. Major antibodies and buffers used in the NADAL® *Trichomonas Vaginalis* Test are as follows. Test line region (T): mouse anti-*Trichomonas* antibodies (26 µg/test); control line region (C): polyclonal goat anti-mouse IgG (62 µg/test); conjugate pad: mouse anti-*Trichomonas* antibodies-colloidal gold conjugate (18 µg/test).

The NADAL® *Trichomonas Vaginalis* Test is a sandwich immunoassay. The test cassette contains a nitrocellulose membrane strip with immobilised mouse anti-*Trichomonas* antibodies in the test line region (T). Goat-anti-mouse antibodies are immobilised in the control line region (C).

The sample is added to the specimen well (S) of the test cassette. The gold conjugate pad contains further mouse anti-*Trichomonas* antibodies which are coupled with colloidal gold. The analyte (*Trichomonas* antigen) in the sample reacts with the colloidal gold coupled mouse anti-*Trichomonas* antibodies of the gold conjugate pad, thus forming antigen-antibody-colloidal gold complexes. Meanwhile the liquid moves along the membrane, transporting these complexes by capillary

action supported by the absorption pad to the respective immobilised mouse anti-*Trichomonas* antibodies. Here they are trapped and form a sandwich complex consisting of immobilised antibody-antigen (analyte)-antibody-colloidal gold. Only when the applied sample contains a certain concentration of *Trichomonas* antigen, will the formation of this sandwich complex result in a visible coloured line in the test line region (T) of the membrane. If no *Trichomonas* antigen is present in the sample, no test line will appear in the test line region of the membrane. The gold conjugates continue to migrate to the control line region (C) of the nitrocellulose membrane. There, these conjugates form complexes with the immobilised goat anti-mouse IgG on the membrane resulting in the formation of a coloured control line (C).

The appearance of a coloured line in the control line region serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 10 NADAL® *Trichomonas Vaginalis* test cassettes (incl. disposable droppers)
- 10 extraction tubes
- Provided additional material according to 93/42/EEC: 10 swabs CE 0197

Dalian Rongbang Medical Healthy Devices Co., Ltd., Guoyi St, Ganjingzi, Dalian, Liaoning, China (authorised EU representative Yao Tong S.L., Ausiàs March, 92; 08013 Barcelona, Spain)

- 1 extraction buffer
- 1 package insert

5. Additional Materials Required

- Timer

6. Storage & Stability

NADAL® *Trichomonas Vaginalis* test kits should be stored at 4-30°C and used until the expiration date indicated on the package. NADAL® *Trichomonas Vaginalis* test cassettes are sensitive to humidity as well as to heat.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Avoid splashing or aerosol formation.
- Clean spills thoroughly using an appropriate disinfectant.

- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).

8. Specimen Collection and Preparation

Specimen Collection

Remove a sterile swab from the pouch. Gently open the vaginal opening and insert the swab about 2 – 3 inches (5 - 7.5 cm) into the vagina. Use extreme care if the patient is pregnant. Improper collection may result in poor visual reading of the test and cause invalid results.

Specimen Preparation

- Bring all samples to room temperature prior to testing.
- Add 7 drops of the extraction buffer to the extraction tube. Put the specimen swab into the tube and swirl it vigorously to mix the reagents for about 15 seconds. Then incubate the mixture with the swab in the tube at room temperature for 10 minutes.
- Swirl the swab vigorously for 15 seconds, then expel as much liquid as possible from the swab by rotating and pressing the fiber tip of the swab against the wall of the tube. Discard the swab. Mix the content of the tube by gently swirling the tube. The swab extract should be tested immediately or stored refrigerated.
- If the swab is not extracted immediately, it can be stored refrigerated, preferably in a transportation tube, for up to 5 days. Do not freeze swabs. Swabs may be transported to the test site under ambient conditions. Transport media should not be used.

9. Test Procedure

Bring tests, buffer and specimens to room temperature (15-30°C) prior to testing. The test cassette should remain in the foil pouch until use.

1. Remove the NADAL® *Trichomonas Vaginalis* test cassette from the foil pouch and place it on a flat, dry surface.
2. Holding a disposable dropper above the sample well (S) of the

test cassette, add 4 drops of the extracted specimen to the sample well (S).

As the test begins to run, you will observe a coloured liquid migrate along the membrane.

3. Read the test result after 15 minutes. Do not interpret test results after more than 15 minutes.

10. Result Interpretation

The above stated result interpretation time is based on reading test results at room temperature (15-30°C). If your room temperature is significantly lower than 15°C, then the interpretation time should be properly increased.

Positive:

Two coloured lines appear in the result window. One line appears in the control line region (C) and the other line appears in the test line region (T).

Note: Generally, the higher the analyte level in the specimen is, the stronger the colour of the “T” line will be. When the analyte level in the specimen is close to but still within the sensitivity limit of the test, the colour of the “T” line will be very faint.

Negative:

Only one coloured line appears in the control line region (C). No coloured line appears in the test line region (T).

Invalid:

The control line fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the invalid result.

Note: A positive result will not change once it has been established after exactly 15 minutes. However, in order to prevent any incorrect results, the test result should not be interpreted after more than 15 minutes.

11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

12. Limitations

- The NADAL® *Trichomonas Vaginalis* Test is limited to the detection of *Trichomonas vaginalis* in swab specimens. Although the test is very accurate in detecting *Trichomonas vaginalis*, a low incidence of false results may occur. Other

clinically available tests are required if questionable results are obtained.

- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

13. Performance Characteristics

Comparison studies were conducted with collected vaginal samples using the NADAL® Trichomonas Vaginalis Test in comparison with traditional microbiology methods and with a commercially available BD Affirm™ Microbial Identification System (USA).

Sensitivity testing of the NADAL® Trichomonas Vaginalis Test in comparison with traditional microbiology methods (n = 300):

NADAL® Trichomonas Vaginalis Test	Traditional microbiology methods			
		+	-	Total
	+	137	4	141
-	13	146	159	
Total	150	150	300	

Relative sensitivity: 91.3% (137/150)

Relative specificity: 97.3% (146/150)

Sensitivity testing of the NADAL® Trichomonas Vaginalis Test in comparison with BD Affirm™ Microbial Identification System (n = 300):

NADAL® Trichomonas Vaginalis Test	BD Affirm™ Microbial Identification System			
		+	-	Total
	+	137	7	144
-	5	151	156	
Total	142	158	300	

Relative sensitivity: 96.5% (137/142)

Relative specificity: 95.6% (151/158)

Cross-reactivity and Interference Studies

Specificity was assessed by studying samples from *Trichomonas* negative subjects, containing no other closely related microorganisms and spiked with similar microorganisms presented in the table below. These samples were tested using the NADAL® Trichomonas Vaginalis Test by a replicate of 10. A sample was classified as negative, when no coloured test line appeared but the coloured control line was visible within 15 minutes. A sample was classified as positive, when both the control and test lines were visible within 15 minutes.

The results are presented in the following table:

Number of tests	Microorganisms spiked with samples	Results of NADAL® Trichomonas Vaginalis Test	
		Spiked negative samples to <i>Trichomonas</i>	Spiked positive samples to <i>Trichomonas</i>
10	<i>Acinetobacter calcoaceticus</i>	All (-)	All (+)
10	<i>Proteus vulgaris</i>	All (-)	All (+)
10	<i>Salmonella typhi</i>	All (-)	All (+)
10	<i>Acinetobacter spp</i>	All (-)	All (+)

Number of tests	Microorganisms spiked with samples	Results of NADAL® Trichomonas Vaginalis Test	
		Spiked negative samples to <i>Trichomonas</i>	Spiked positive samples to <i>Trichomonas</i>
10	<i>Staphylococcus aureus</i>	All (-)	All (+)
10	<i>Escherichia coli</i>	All (-)	All (+)
10	<i>Moraxella catarrhalis</i>	All (-)	All (+)
10	<i>Gardnerella vaginalis</i>	All (-)	All (+)
10	<i>Streptococcus faecalis</i>	All (-)	All (+)
10	<i>Neisseria lactamica</i>	All (-)	All (+)
10	<i>Streptococcus faecium</i>	All (-)	All (+)
10	<i>Pseudomonas aeruginosa</i>	All (-)	All (+)
10	<i>Candida albicans</i>	All (-)	All (+)

Potentially interfering chemical substances (presented in the table below) were supplemented to clinically defined negative normal specimens. The samples were tested using NADAL® Trichomonas Vaginalis Test by a replicate of 5. A sample was classified as negative, when no coloured test line appeared but the coloured control line was visible within 15 minutes. A sample was classified as positive, when both the control and test lines were visible within 15 minutes.

Number of tests done	Potentially interfering chemical substances spiked with negative samples	
5	Acetaminophen, 20 mg/dl	Acetyl salicylic acid, 20 mg/dl
5	Ascorbic acid, 20 mg/dl	Atropine, 20 mg/dl
5	Bilirubin, 60 mg/dl	Caffeine, 20 mg/dl
5	Creatinine, 20 mg/dl	Gentamic acid, 20 mg/dl
5	Glucose, 2,000 mg/dl	Haemoglobin, 500 mg/dl
5	Ketones, 40 mg/dl	Mestranol, 3 mg/dl
5	Nitrite, 20 mg/dl	Penicillin, 40,000 U/dl
5	Sodium heparin, 3 mg/dl	Lithium heparin, 3 mg/dl

A set of interfering chemical substances (presented in the table below) was selected for a high dose interference study. Clinically defined *Trichomonas* negative samples were spiked with the set of selected high dose interference substances and tested using the NADAL® Trichomonas Vaginalis Test.

High dose interference study		
Interference Substance	Concentration	Result (3 lots)
Triglyceride	500 mg/ml	-/-
Bilirubin	0.1 mg/ml	-/-
Prostatic acid phosphatase	1000 mIU	-/-
Albumin	20 mg/ml	-/-

None of the tested substances (presented in the table above) showed any interference with clinically defined negative specimens. The negative samples spiked with potentially

interfering substances showed consistently negative test results.

Reproducibility Study

The reproducibility study was carried out on 3 different sites. Negative, weak positive and strong positive samples were tested using 120 NADAL® Trichomonas Vaginalis test kits in total. The samples were tested twice a day, in 2 different assays, daily for 20 days. This allowed for separate tests of between-day, between-assay and within-day, which showed consistent results.

Reproducibility study based on study sites

Trichomonas test readings (number of tests)			
Test sites	Negative	Weak positive	Strong positive
1	N (40)	P (40)	P (40)
2	N (40)	P (40)	P (40)
3	N (40)	P (40)	P (40)

The test results on the 3 sites were in 100% agreement.

14. References

1. Bowden FJ, Garnett GP (2000) Trichomonas vaginalis epidemiology: parameterising and analysis, a model of treatment intervention. Sex Trans Infect 76:248-256.
2. Wang CC, McClellandRS, ReillyM, OverbaughJ, EmerySR, MandaliyaK, Chohan B, Ndinya Acholaj, Bwayoj,Kreiss JK (2001) The effect of treatment of vaginal infection on shedding of immunodeficiency virus type 1. J Infect Dis 183:1017-1022.
3. World Health Organisation (2001) Global prevalence and incidence of selected curable sexually transmitted infections: overview and estimates. World Health Organisation, Geneva, Switzerland.
4. Ledru S, Meda N, Fofana, M, Soula G, Bazie AJ, Chiron JP (1992) Etiologic study of genitourinary infections in women of child bearing age in Bobo-Dioulasso, Burkina FasoSex Transm Dis; 23(2):151-156.
5. Yereli K, Balcioglu IC, Degerli K, Ozbilgin A, Daldal N (1997) Incidence of Trichomonas vaginalis among women having vaginal discharge, in Manisa, Turkey J Egypt Soc Parasitol 27 (3): 905 – 911.

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14. Bibliographie

- Bowden FJ, Garnett GP (2000) Trichomonas vaginalis epidemiology: parameterising and analysis, a model of treatment intervention. Sex Trans Infect 76:248-256.
- Wang CC, McClellandRS, ReillyM, OverbaughJ, EmerySR, MandalayK, Chohan B, Ndinya Acholaj, BwayoJ, Kreiss JK (2001) The effect of treatment of vaginal infection on shedding of immunodeficiency virus type 1. J Infect Dis 183:1017-1022.
- World Health Organisation (2001) Global prevalence and incidence of selected curable sexually transmitted infections: overview and estimates. World Health Organisation, Geneva, Switzerland.

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14. Referencias

- Bowden FJ, Garnett GP (2000) Trichomonas vaginalis epidemiology: parameterising and analysis, a model of treatment intervention. Sex Trans Infect 76:248-256.
- Wang CC, McClellandRS, ReillyM, OverbaughJ, EmerySR, MandaliyaK, Chohan B, Ndinya Aholaj, Bwayoj,Kreiss JK (2001) The effect of treatment of vaginal infection on shedding of immunodeficiency virus type 1. J Infect Dis 183:1017-1022.
- World Health Organisation (2001) Global prevalence and incidence of selected curable sexually transmitted infections: overview and estimates. World Health Organisation, Geneva, Switzerland.
- Ledru S, Meda N, Fofana, M, Soula G, Bazie AJ, Chiron JP (1992) Etiologic study of genitourinary infections in women of child bearing age in Bobo-Dioulasso, Burkina FasoSex Transm Dis; 23(2):151-156.
- Yereli K, Balcioglu IC, Degerli K, Ozbilgin A, Dalal N (1997) Incidence of Trichomonas vaginalis among women having vaginal discharge, in Manisa, Turkey J Egypt Soc Parasitol 27 (3): 905 – 911.

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- World Health Organisation (2001) Global prevalence and incidence of selected curable sexually transmitted infections: overview and estimates. World Health Organisation, Geneva, Switzerland.
- Ledru S, Meda N, Fofana, M, Soula G, Bazie AJ, Chiron JP (1992) Etiologic study of genitourinary infections in women of child bearing age in Bobo-Dioulasso, Burkina Faso Sex Transm Dis; 23(2):151-156.
- Yereli K, Balcioglu IC, Degerli K, Ozbilgin A, Daldal N (1997) Incidence of Trichomonas vaginalis among women having vaginal discharge, in Manisa, Turkey J Egypt Soc Parasitol 27 (3): 905 – 911.

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14. Reference

- Bowden FJ, Garnett GP (2000) Trichomonas vaginalis epidemiology: parameterising and analysis, a model of treatment intervention. Sex Trans Infect 76:248-256.
- Wang CC, McClellandRS, ReillyM, OverbaughJ, EmerySR, MandaliyaK, Chohan B, Ndinya Acholaj, Bwayoj,Kreiss JK (2001) The effect of treatment of vaginal infection on shedding of immunodeficiency virus type 1. J Infect Dis 183:1017-1022.
- World Health Organisation (2001) Global prevalence and incidence of selected curable sexually transmitted infections: overview and estimates. World Health Organisation, Geneva, Switzerland.
- Ledru S, Meda N, Fofana, M, Soula G, Bazie AJ, Chiron JP (1992) Etiologic study of genitourinary infections in women of child bearing age in Bobo-Dioulasso, Burkina FasoSex Transm Dis; 23(2):151-156.
- Yereli K, Balcioglu IC, Degerli K, Ozbilgin A, Daldal N (1997) Incidence of Trichomonas vaginalis among women having vaginal discharge, in Manisa, Turkey J Egypt Soc Parasitol 27 (3): 905 – 911.

Rev. 0, 2015-08-11 TF

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conforme aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	<i>in-vitro</i> -Diagnostika	<i>in-vitro</i> diagnostic medical device	Dispositif médical de diagnostic <i>in-vitro</i>	Producto sanitario para diagnóstico <i>in-vitro</i>	Dispositivo medico-diagnostico <i>in-vitro</i>	Tylko do diagnostyki <i>in-vitro</i>
	Temperaturbegrenzung	Temperature limitation	Limites de température	Limite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Numéro de lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico <i>in-vitro</i>	Diagnostický zdravotnický prostředek <i>in-vitro</i>	<i>in-vitro</i> - diagnostiikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för <i>in-vitro</i> -diagnostik	Medisch hulpmiddel voor <i>in-vitro</i> diagnostiek	Medicinsk udstyr til <i>in-vitro</i> -diagnostik	<i>in-vitro</i> diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperaturbegränsning	Temperatuurlimiet	Temperaturbegrænsning	Temperaturbegrensning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merking
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebuje do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Best il l ingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrekkelig for <n> tester

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0

Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0

Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29

Free Tel: 0800 291 565

Fax: +49 290 10-50

Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18

Free Tel – UK: 0808 234 1237

Free Tel – IRE: 1800 555 080

Fax: +49 290 10-50

France:

France Tel: 0800 915 240

France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720

Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82

Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16

Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759

Free Tel: 900 938 315

Fax: +49 941 290 10-50

Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34

Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44

Free Tel: 00 800 491 15 95

Fax: +49 941 290 10-50

Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735

Tel. Verde: 800 849 230

Fax: +49 941 290 10-50

Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600

Free Tel: 0800 0222 890

Fax: +31 70 30 30 775

Free Fax: 0800 024 9519

Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607

Free Tel: +45 80 88 87 53

Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40

Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH