

TRACEABLE VALUE FOR STANDARDS

| COMPONENT | Cat. Number | TRACEABILITY | Matrix |
|-----------------------------------|---------------|---|---------------------------|
| CHOLESTEROL HDL/LDL DIRECT | 11693 | CDC Reference Measurement Procedure (Center for Disease Control and Prevention) | Human serum |
| BILIRUBIN | 11513 | Standard Reference Material 916a (National Institute of Standards and Technology, USA) | Aqueous |
| ADENOSINE DEAMINASE (ADA) | 18052 | Reference Material BCR-647 (IRMM) | Bovine ADA |
| CARBON DIOXIDE (CO ₂) | 11833 / 11822 | Standard Reference Material SRM 351a (National Institute of Standards and Technology, NIST) | Sodium Hydrogen Carbonate |
| PROTEIN(URINE) STANDARD | 11603 | Standard Reference Material SRM 927 (National Institute of Standards and Technology, NIST) | Bovine Albumin |



- Certified Management System
- EN ISO 9001
- EN ISO 13485